



Rapid Testing. Reliable Results.

Legionella Sampling Packing Instructions

Small Box (White Kraft 8W x 8D x 10.5H)

1 Sealable Insulated Bag

Option 1: Maximum of 18–250 mL Bottles

Option 2: Maximum of 4–1L Bottles



Medium Box (Plain Kraft 16W x 16D x 9H)

1 Sealable Insulated Bag

Option 1: Maximum of 72–250 mL Bottles (18 per bag)

Option 2: Maximum of 16–1L Bottles (4 per bag)



Large Bulk à la carte (28W x 18D x 18H)

of Small Boxes (White Kraft) - flat packed

of Medium Kit Boxes (Plain Kraft) - flat packed

of Sealable Insulated Bags - flat packed

of 250 mL Bottles

of 1L Bottles

Packing Instructions (Please Follow Carefully)

1. Insert sealable insulated bag into box
2. Load bag with bottles (250 mL and/or 1 L)
3. Remove cellophane and seal insulated bag
4. If headspace exists, fill the gap with adequate packing materials for secure shipment
5. Print off a copy of Chain of Custody/Sample Submission Form and place in package
6. Securely tape box closed for shipment
7. Attach shipping label (UPS/FedEx)
8. Transport box to carrier location (UPS/FedEx)



Rapid Testing. Reliable Results.

Sampling Guidance Disclaimer. *The sampling guidance and tips provided by Q Laboratories are intended for general informational purposes only and do not constitute project-specific instructions, professional advice, or a substitute for site-specific procedures or regulatory requirements. Sampling conditions, objectives, and regulatory obligations vary, and this guidance may not be applicable to all situations. Users are solely responsible for determining the suitability of any sampling method for their specific application and for ensuring compliance with all applicable laws, regulations, standards, and contractual obligations. Q Laboratories makes no representations or warranties, express or implied, regarding the accuracy or completeness of this information and assumes no responsibility or liability for outcomes resulting from its use or reliance. For project-specific guidance, users should consult applicable regulations or qualified professionals.*